

## T9982 MPGEL Re-enterable Clear Sealant 1 litre

MPGEL Re-enterable Clear Sealant which is mixed in the jug provided for pouring into ground-burial joints and enclosures requiring sealing against wicking. Sealant can be removed in case of future maintenance requirements.



### TECHNICAL SPECIFICATION

- Curing time @23c : 20 minutes
- No halogens and plasticizers
- Low curing temperature
- Mixing time 8 minutes
- 1ltr (2 x 500ml 50:50 mix)

### Safety Notes:

- Personal protective clothing (gloves, eye protection) should be worn
- Ensure adequate ventilation is available
- Mild irritation may occur should contact with skin happen, in which case immediately wash affected area with soap and water, seeking medical advice if reaction persists
- If in contact with eyes, bathe the eye with running water for 15 minutes, and seek medical attention if discomfort continues

### COMPLIANCE & WARRANTY

- 2 Year Warranty
- CE / RoHS Compliant